



## ABCEPTA CUSTOM SERVICES:

### Transfection Protocol

#### Materials

- u Cell line cultured in the appropriate growth medium
- u Plasmid DNA
- u Lipofectamine
- u Basic Medium (without serum)
- u Appropriate cell culture plates and supplies

#### Procedure

1. Adherent cells: Plate  $0.5-2 \times 10^5$  cells in 500 uL of growth medium without antibiotics at the day before transfection, so that cells will be 90-95% confluent at the time of transfection.  
Suspension cells: Just plate  $4-8 \times 10^5$  cells in 500 uL of growth medium without antibiotics.
2. For each transfection sample, prepare complexes as follows:
  - a. Dilute DNA in 50 uL of Basic Medium. Mix gently.
  - b. Mix Lipofectamine gently before use, then dilute the appropriate amount (DNA: Lipofect = 1:2.5) in 50 uL of Basic Medium. Incubate for 5 minutes at room temperature.
  - c. After the 5 minute incubation, combine the diluted DNA with diluted Lipofectamine (total volume = 100 uL). Mix gently and incubate for 20 minutes at room temperature (solution may appear cloudy).  
**Note: All the steps should be finished with in 25 minutes.**
3. Add the 100 uL of complexes to each well containing cells and medium. Mix gently by shaking the plate left and right.
4. Incubate cells at 37 °C in a CO<sub>2</sub> incubator for 18-48 hours before testing for protein expression. Medium may be changed after 4-6 hours.
5. Stable cell lines: Passage cells at a 1:10 (or higher dilution) into fresh growth medium 24 hours after transfection. Add selective medium (if desired) the following day.