

# **Biotinylated Anti-Human ApoE3 Antibody**

Catalog # ABG10011

# **Specification**

# Biotinylated Anti-Human ApoE3 Antibody - Product Information

Application WB, E
Reactivity Human
Host Rabbit
Clonality Polyclonal

## Biotinylated Anti-Human ApoE3 Antibody - Additional Information

## **Preparation**

Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant hApoE3. Anti-Human ApoE3 specific antibody was purified by affinity chromatography and then biotinylated. <strong>\*Cross Reactivity: </strong> human ApoE2, human ApoE4 in both wester

#### WesternBlot

To detect hApoE3 by Western Blot analysis this antibody can be used at a concentration of  $0.1 - 0.2 \,\mu g/ml$ . Used in conjunction with compatible secondary reagents the detection limit for recombinant hApoE3 is  $1.5 - 3.0 \, ng/lane$ , under either reducing or non-reducing

## Sandwich

To detect hApoE3 by sandwich ELISA (using 100  $\mu$ l/well antibody solution) a concentration of 0.25 – 1.0  $\mu$ g/ml of this antibody is required. This biotinylated polyclonal antibody, in conjunction with BioGems' Polyclonal Anti-Human ApoE3 (60-062P) as a capture antibody, allows the detection of at least 0.2 – 0.4 ng/well of recombinant hApoE3.

## **Direct**

To detect hApoE3 by direct ELISA (using 100  $\mu$ l/well antibody solution) a concentration of 0.25 – 1.0  $\mu$ g/ml of this antibody is required. This biotinylated polyclonal antibody, in conjunction with compatible secondary reagents, allows the detection of at least 0.2 – 0.4 ng/well of recombinant hApoE3.

#### **Formulation**

A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.

#### Reconstitution

Centrifuge vial prior to opening. Reconstitute in sterile PBS containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml.

## **Storage**

-20°C

## **Precautions**

Biotinylated Anti-Human ApoE3 Antibody is for research use only and not for use in diagnostic or therapeutic procedures.



Tel: 858.875.1900 Fax: 858.875.1999



# **Biotinylated Anti-Human ApoE3 Antibody - Protocols**

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- <u>Immunofluorescence</u>
- <u>Immunoprecipitation</u>
- Flow Cytomety
- Cell Culture

**Biotinylated Anti-Human ApoE3 Antibody - Images**