

Anti-Murine B7-1 Antibody

Catalog # ABG10016

Specification

Anti-Murine B7-1 Antibody - Product Information

Application WB, E
Reactivity Mouse
Host Rabbit
Clonality Polyclonal

Anti-Murine B7-1 Antibody - Additional Information

Preparation

Produced from sera of rabbits immunized with highly pure recombinant Murine B7-1. Anti-Murine B7-1 specific antibody was purified by affinity chromatography employing an immobilized Murine B7-1 matrix.

WesternBlot

To detect mB7-1 by Western Blot analysis this antibody can be used at a concentration of 0.1- 0.2 μ g/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant mB7-1 is 1.5-3.0 ng/lane, under either reducing or non-reducing conditions.

Sandwich

To detect mB7-1 by sandwich ELISA (using 100 μ l/well antibody solution) a concentration of 0.5 - 2.0 μ g/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with BioGems' Biotinylated Anti-Murine B7-1 (61-401BT) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant mB7-1.

Formulation

A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.

Reconstitution

Centrifuge vial prior to opening. Reconstitute in sterile water to a concentration of 0.1-1.0 mg/ml.

Storage

-20°C

Precautions

Anti-Murine B7-1 Antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Anti-Murine B7-1 Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot





• <u>Immunohistochemistry</u>

- <u>Immunofluorescence</u>
- <u>Immunoprecipitation</u>
- Flow Cytomety
- Cell Culture

Anti-Murine B7-1 Antibody - Images