

Anti-Human Eotaxin-2 Antibody

Catalog # ABG10091

Specification

Anti-Human Eotaxin-2 Antibody - Product Information

Application WB, IHC, E
Reactivity Human
Host Mouse
Clonality Monoclonal

Anti-Human Eotaxin-2 Antibody - Additional Information

Preparation

Produced in BALB/c mice using highly pure (>98%) recombinant human Eotaxin-2 as the immunizing antigen. This IgG1 < sub > K < / sub > antibody was purified from ascites fluid by antigen affinity chromatography.

WesternBlot

To detect hEotaxin-2 by Western Blot analysis this antibody can be used at a concentration of 0.25- $0.50 \mu g/ml$. Used in conjunction with compatible secondary reagents the detection limit for recombinant hEotaxin-2 is 0.5- $1.0 \mu g/lane$, under non-reducing conditions.

Sandwich

In a sandwich ELISA (assuming 100μ I/well), a concentration of 2.0-4.0 µg/ml of this antibody will detect at least 80.0 pg/well of recombinant human Eotaxin-2 when used with BioGems' biotinylated antigen affinity purified anti-human Eotaxin-2 (60-123BT) as the detection antibody at a concentration of approximately 1.0-2.0 µg/ml.

Immunohistochemistry

This antibody stained formalin-fixed, paraffin-embedded sections of human tonsil. The recommended concentration is $10.0~\mu g/ml$ with an overnight incubation at $4^{\circ}C$. An HRP-labeled polymer detection system was used with a DAB chromogen. Heat induced antigen retrieval with a pH 8.0~EDTA buffer is recommended. Optimal concentrations and conditions may vary. Protocol and staining provided by Dr. Lauren Binge, Laboratory of Prof. Charles Mackay, Monash University, Australia.

Formulation

A sterile filtered antibody solution was lyophilized from PBS.

Reconstitution

Centrifuge vial prior to opening. Reconstitute in sterile water to a concentration of 0.1-1.0 mg/ml.

Storage

-20°C

Precautions

Anti-Human Eotaxin-2 Antibody is for research use only and not for use in diagnostic or therapeutic procedures.



Tel: 858.875.1900 Fax: 858.875.1999



Anti-Human Eotaxin-2 Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- <u>Immunofluorescence</u>
- <u>Immunoprecipitation</u>
- Flow Cytomety
- Cell Culture

Anti-Human Eotaxin-2 Antibody - Images