

Anti-Human GRO-y Antibody

Catalog # ABG10151

Specification

Anti-Human GRO-y Antibody - Product Information

Application WB, IHC, E
Reactivity Human
Host Rabbit
Clonality Polyclonal

Anti-Human GRO-y Antibody - Additional Information

Preparation

Produced from sera of rabbits pre-immunized with highly pure recombinant Human GRO- γ . Anti-Human GRO- γ specific antibody was purified by affinity chromatography employing immobilized Human GRO- γ matrix.

WesternBlot

To detect Human GRO- γ by Western Blot analysis this antibody can be used at a concentration of 0.1-0.2 μ g/ml. When used in conjunction with compatible secondary reagents, the detection limit for recombinant Human GRO- γ is 1.5-3.0 ng/lane, under either reducing or non-reducing conditions.

Sandwich

To detect Human GRO- γ by sandwich ELISA (using 100 μ l/well antibody solution) a concentration of 0.5 - 2.0 μ g/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with BioGems' Biotinylated Anti-Human GRO- γ (60-172BT) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant Human GRO- γ .

Immunohistochemistry

This antibody stained formalin-fixed, paraffin-embedded sections of human colon/rectum adenocarcinoma. The recommended concentration is 0.7 mg/ml-0.5mg/ml with an overnight incubation at 4°C. An HRP-labeled polymer detection system was used with a DAB chromogen. Heat induced antigen retrieval with a pH 6.0 sodium citrate buffer is recommended. Optimal concentrations and conditions may vary.

Neutralization

To yield one-half maximal inhibition [ND_{>50}] of the biological activity of Human GRO- γ (100 ng/ml), a concentration of 1.0-2.0 µg/ml of this antibody is required.

Formulation

A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.

Reconstitution

Centrifuge vial prior to opening. Reconstitute in sterile water to a concentration of 0.1-1.0 mg/ml.

Storage

-20°C



Tel: 858.875.1900 Fax: 858.875.1999

Precautions

Anti-Human GRO-y Antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Anti-Human GRO-y Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- <u>Immunofluorescence</u>
- Immunoprecipitation
- Flow Cytomety
- Cell Culture

Anti-Human GRO-γ Antibody - Images