

**His-Tag Antibody**  
**Rabbit Polyclonal Antibody**  
**Catalog # ABV10768****Specification**

---

**His-Tag Antibody - Product Information**

Application	<b>WB, ICC, IP</b>
Reactivity	<b>All Species</b>
Host	<b>Rabbit</b>
Clonality	<b>Polyclonal</b>
Isotype	<b>Rabbit IgG</b>

**His-Tag Antibody - Additional Information**

Application & Usage	<b>Western blotting (1 µg/ml), immunoprecipitation (40 µg/ml), and immunocytochemistry (2 µg/ml). However, the optimal conditions should be determined individually. The antibody detects His-tagged recombinant proteins.</b>
---------------------	--

**Other Names**

His tag, HHHHHH tag

**Target/Specificity**

His-Tag

**Antibody Form**

Liquid

**Appearance**

Colorless liquid

**Formulation**

100 µg (0.5 mg/ml) peptide affinity purified rabbit anti-His-Tag polyclonal antibody in phosphate buffered saline (PBS, pH 7.2), containing 50% glycerol, 1% BSA, 0.02% thimerosal.

**Handling**

The antibody solution should be gently mixed before use.

**Reconstitution & Storage**

-20 °C

**Background Descriptions****Precautions**

His-Tag Antibody is for research use only and not for use in diagnostic or therapeutic procedures.

**His-Tag Antibody - Protein Information**

## **His-Tag Antibody - Protocols**

Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)

## **His-Tag Antibody - Images**

## **His-Tag Antibody - Background**

A variety of plasmids contain DNA that encodes an N-terminal tag consisting of six histidine (His) residues, followed by an extended multiple cloning sites. The 6xHis tag on the expressed recombinant proteins allows for efficient coupling to Ni<sup>++</sup> affinity resins and purification by a single step chromatography. As is the case with other protein tags, this polyhistidine tag can be cleaved at sites recognized by enzymes such as thrombin and enterokinases to isolate the protein of interest.