

## CLASP1 Rabbit mAb

**Catalog # AP75268** 

## Specification

#### **CLASP1** Rabbit mAb - Product Information

**Application Primary Accession** Reactivity Host Clonality

Calculated MW

WB, IHC-P, IHC-F, ICC

**Q7Z460** Human Rabbit

**Monoclonal Antibody** 

169451

#### CLASP1 Rabbit mAb - Additional Information

Gene ID 23332

**Other Names** CLASP1

**Dilution** 

WB~~1/500-1/1000 IHC-P~~N/A IHC-F~~N/A

ICC~~N/A

**Format** Liquid

### **CLASP1** Rabbit mAb - Protein Information

Name CLASP1

Synonyms KIAA0622, MAST1

#### **Function**

Microtubule plus-end tracking protein that promotes the stabilization of dynamic microtubules. Involved in the nucleation of noncentrosomal microtubules originating from the trans-Golgi network (TGN). Required for the polarization of the cytoplasmic microtubule arrays in migrating cells towards the leading edge of the cell. May act at the cell cortex to enhance the frequency of rescue of depolymerizing microtubules by attaching their plus-ends to cortical platforms composed of ERC1 and PHLDB2. This cortical microtubule stabilizing activity is regulated at least in part by phosphatidylinositol 3-kinase signaling. Also performs a similar stabilizing function at the kinetochore which is essential for the bipolar alignment of chromosomes on the mitotic spindle.

## **Cellular Location**

Cytoplasm, cytoskeleton. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome. Chromosome, centromere, kinetochore Cytoplasm, cytoskeleton, spindle. Golgi apparatus, trans-Golgi network. Note=Localizes to microtubule plus ends. Localizes to centrosomes,





Tel: 858.875.1900 Fax: 858.875.1999

kinetochores and the mitotic spindle from prometaphase Subsequently localizes to the spindle midzone from anaphase and to the midbody from telophase. In migrating cells localizes to the plus ends of microtubules within the cell body and to the entire microtubule lattice within the lamella. Localizes to the cell cortex and this requires ERC1 and PHLDB2

### **CLASP1** Rabbit mAb - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- <u>Immunofluorescence</u>
- <u>Immunoprecipitation</u>
- Flow Cytomety
- Cell Culture

# CLASP1 Rabbit mAb - Images

