

### **Chicken IgG Texas Red™**

Catalog # ASR1262

#### **Specification**

Conjugate

# **Chicken IgG Texas Red™ - Product Information**

Description CHICKEN IgG whole molecule Texas Red™

conjugated Texas Red®

FP Value 3.0 moles Texas Red® per mole of Chicken

IaG

Physical State Lyophilized

Host Isotype IgG

Buffer 0.02 M Potassium Phosphate, 0.15 M

**Sodium Chloride, pH 7.2** 

Species of Origin
Reconstitution Volume
Chicken
1.0 mL

Reconstitution Buffer Restore with deionized water (or

equivalent)

Stabilizer 10 mg/mL Bovine Serum Albumin (BSA) -

Immunoglobulin and Protease free

Preservative 0.01% (w/v) Sodium Azide

# **Chicken IgG Texas Red™ - Additional Information**

#### **Shipping Condition**

**Ambient** 

#### **Purity**

This product was prepared from normal serum by delipidation, salt fractionation, ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Chicken IgG and anti-Chicken Serum.

#### **Storage Condition**

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

# **Precautions Note**

This product is for research use only and is not intended for therapeutic or diagnostic applications.

#### **Chicken IgG Texas Red™ - Protein Information**


#### Chicken IgG Texas Red™ - Protocols



Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- Cell Culture

# Chicken IgG Texas Red™ - Images



SDS-PAGE of Chicken IgG Whole Molecule Texas Red  $^{\mathsf{m}}$  Conjugated . Lane M: 5  $\mu$ L Opal Prestained Marker . Lane 1: Reduced Chicken IgG Whole Molecule Texas Red  $^{\mathsf{m}}$  Conjugated . Lane 2: Reduced Chicken IgG F(c) Fragment . Lane 3: Reduced Chicken IgG F(ab) Fragment . Lane 4: Reduced Chicken IgM Whole Molecule . Load: 1  $\mu$ g per lane. Predicted/Observed size: IgG at 72 and 25 kDa; F(c) at 25 kDa; F(ab) at 25 kDa; IgM at 75 kDa. Observed F(c) Fragment migrates slightly higher. Other bands: Chicken IgG heavy chain alternative splicing variant at approximately 40 kDa in Lane 1.

# Chicken IgG Texas Red™ - Background

This product is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.