

Rat IgG Fab Peroxidase

Catalog # ASR1746

Specification

Rat IgG Fab Peroxidase - Product Information

Description RAT IgG F(ab) fragment Peroxidase

conjugated

Conjugate Peroxidase (Horseradish)

Physical State
Host Isotype
Buffer
Lyophilized
IgG F(ab)
0.01 M Sod

0.01 M Sodium Phosphate, 0.15 M Sodium

Chloride, pH 7.2

Species of Origin
Reconstitution Volume
Rat
1.0 mL

Reconstitution Buffer Restore with deionized water (or

equivalent)

Stabilizer 10 mg/mL Bovine Serum Albumin (BSA) -

Immunoglobulin and Protease free

Preservative 0.01% (w/v) Thimerosal

Rat IgG Fab Peroxidase - Additional Information

Shipping Condition

Ambient

Purity

This product was prepared from normal serum by delipidation, salt fractionation, ion exchange chromatography followed by papain digestion and extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti- Rat $\lg G$, anti-Rat $\lg G$ F(ab')2 and anti-Rat Serum. No reaction was observed against anti-Rat $\lg G$ F(c) or anti-Papain.

Storage Condition

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Precautions Note

This product is for research use only and is not intended for therapeutic or diagnostic applications.

Rat IgG Fab Peroxidase - Protein Information

Rat IgG Fab Peroxidase - Protocols

Provided below are standard protocols that you may find useful for product applications.






Tel: 858.875.1900 Fax: 858.875.1999

- Western Blot
- Blocking Peptides
- Dot Blot
- Immunohistochemistry
- Immunofluorescence
- <u>Immunoprecipitation</u>
- Flow Cytomety
- Cell Culture

Rat IgG Fab Peroxidase - Images



SDS-PAGE of Rat IgG F(ab) Fragment Peroxidase Conjugated . Lane M: 3 µL Opal Prestained Marker . Lane 1: Reduced Rat IgG Whole Molecule . Lane 2: Reduced Rat IgG F(c) Fragment . Lane 3: Reduced Rat IgG F(ab) Fragment Peroxidase Conjugated . Lane 4: Reduced Rat IgM Whole Molecule . Load: 1 μg of IgG, F(c), F(ab); 1.5 μg of IgM. Predicted/Observed size: IgG at 55 and 25 kDa; F(c) at 25 kDa; F(ab) at 25 kDa; IgM at 78 and 25 kDa. Observed F(c) Fragment migrates slightly higher.