

Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody
Goat Polyclonal, Peroxidase (Horseradish)
Catalog # ASR1884**Specification****Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Product Information**

| | |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------|
| Description | Anti-CHICKEN IgG (H&L) (GOAT) Antibody Peroxidase Conjugated |
| Host | Goat |
| Conjugate | Peroxidase (Horseradish) |
| Target Species | Chicken |
| Reactivity | Chicken |
| Clonality | Polyclonal |
| Application | ,1,10,15, |
| Application Note | ELISA 1:10,000-1:150,000;Western Blot 1:1,000-1:5,000;Immunochemistry 1:1,000-1:10,000 |
| Physical State | Lyophilized |
| Host Isotype | IgG |
| Target Isotype | IgG (H&L) |
| Buffer | 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 |
| Immunogen | Anti-Chicken IgG (H&L) was produced by repeated immunization with Chicken IgG whole molecule in goat. |
| Reconstitution Volume | 1.0 mL |
| Reconstitution Buffer | Restore with deionized water (or equivalent) |
| Stabilizer | 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free |
| Preservative | 0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide! |

Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Additional Information**Shipping Condition**

Ambient

Purity

This product was prepared from monospecific antiserum by immunoaffinity chromatography using Chicken IgG coupled to agarose beads. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Goat Serum, Chicken IgG and Chicken Serum.

Storage Condition

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Precautions Note

This product is for research use only and is not intended for therapeutic or diagnostic applications.

Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Protein Information**Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Protocols**

Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)

Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Images**Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Background**

Anti-Chicken IgG antibody generated in goat detects specifically chicken IgG. This secondary antibody anti-Chicken is ideal for investigators who routinely perform ELISA, Sandwich ELISA, titration assays, western-blot, immunoprecipitation and more generally immunoassays.