

### Anti-Horse IgG (H&L) (Alkaline Phosphatase Conjugated) Secondary Antibody

Goat Polyclonal, Alkaline Phosphatase (Calf Intestine)
Catalog # ASR2334

#### **Specification**

# Anti-Horse IgG (H&L) (Alkaline Phosphatase Conjugated) Secondary Antibody - Product Information

Description Anti-HORSE IgG (H&L) (GOAT) Antibody

**Alkaline Phosphatase Conjugated** 

Host Goat

Conjugate Alkaline Phosphatase (Calf Intestine)

Target Species Horse
Clonality Polyclonal
Application WB, IHC, E

Application Note ELISA 1:2,000-1:10,000; Western Blot

1:500-1:2,500;Immunohistochemistry

1:200-1:1,000

Physical State Liquid (sterile filtered)

Host Isotype IgG

Target Isotype IgG (H&L)

Buffer 0.05 M Tris Chloride, 0.15M Sodium

Chloride, 0.001M Magnesium Chloride, 0.0001M Zinc Chloride, 50% (v/v) Glycerol;

pH 8.0

Immunogen Horse IgG whole molecule

Stabilizer 10 mg/mL Bovine Serum Albumin (BSA) -

Immunoglobulin and Protease free

Preservative 0.01% (w/v) Sodium Azide

# Anti-Horse IgG (H&L) (Alkaline Phosphatase Conjugated) Secondary Antibody - Additional Information

#### **Shipping Condition**

Wet Ice

### **Purity**

This product was prepared from monospecific antiserum by immunoaffinity chromatography using Horse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Alkaline Phosphatase (calf intestine), anti-Goat Serum, Horse IgG and Horse Serum.

### **Storage Condition**

Store vial at 4° C before opening. DO NOT FREEZE. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity.

#### **Precautions Note**

This product is for research use only and is not intended for therapeutic or diagnostic applications.



## Anti-Horse IgG (H&L) (Alkaline Phosphatase Conjugated) Secondary Antibody - Protein Information

## Anti-Horse IgG (H&L) (Alkaline Phosphatase Conjugated) Secondary Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- Cell Culture

Anti-Horse IgG (H&L) (Alkaline Phosphatase Conjugated) Secondary Antibody - Images