

Anti-Horse IgG (H&L) (Alkaline Phosphatase Conjugated) Secondary Antibody

Rabbit Polyclonal, Alkaline Phosphatase (Calf Intestine)
Catalog # ASR2460

Specification

Anti-Horse IgG (H&L) (Alkaline Phosphatase Conjugated) Secondary Antibody - Product Information

Description Anti-HORSE IgG (H&L) (RABBIT) Antibody

Alkaline Phosphatase Conjugated

Host Rabbit

Conjugate Alkaline Phosphatase (Calf Intestine)

Target Species Horse
Clonality Polyclonal
Application ,1,10,15,

Application Note ELISA 1:2,000-1:20,000; Western Blot

1:500-1:2,500; Immunochemistry

1:200-1:1,000

Physical State Liquid (sterile filtered)

Host Isotype IgG

Target Isotype IgG (H&L)

Buffer 0.05 M Tris Chloride, 0.15M Sodium

Chloride, 0.001M Magnesium Chloride, 0.0001M Zinc Chloride, 50% (v/v) Glycerol;

pH 8.0

Immunogen Horse IgG whole molecule

Stabilizer 10 mg/mL Bovine Serum Albumin (BSA) -

Immunoglobulin and Protease free

Preservative 0.1% (w/v) Sodium Azide

Anti-Horse IgG (H&L) (Alkaline Phosphatase Conjugated) Secondary Antibody - Additional Information

Shipping Condition

Wet Ice

Purity

This product was prepared from monospecific antiserum by immunoaffinity chromatography using Horse IgG coupled to agarose beads. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Alkaline Phosphatase (calf intestine), anti-Rabbit Serum, Horse IgG and Horse Serum.

Storage Condition

Store vial at 4° C before opening. DO NOT FREEZE. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity.

Precautions Note

This product is for research use only and is not intended for therapeutic or diagnostic applications.



Anti-Horse IgG (H&L) (Alkaline Phosphatase Conjugated) Secondary Antibody - Protein Information

Anti-Horse IgG (H&L) (Alkaline Phosphatase Conjugated) Secondary Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- Cell Culture

Anti-Horse IgG (H&L) (Alkaline Phosphatase Conjugated) Secondary Antibody - Images