

RGN Antibody (Center) Blocking Peptide

Synthetic peptide Catalog # BP7487c

Specification

RGN Antibody (Center) Blocking Peptide - Product Information

Primary Accession

Q15493

RGN Antibody (Center) Blocking Peptide - Additional Information

Gene ID 9104

Other Names

Regucalcin, RC, Gluconolactonase, GNL, Senescence marker protein 30, SMP-30, RGN, SMP30

Target/Specificity

The synthetic peptide sequence used to generate the antibody AP7487c was selected from the Center region of human RGN. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.

Format

Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.

Storage

Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.

Precautions

This product is for research use only. Not for use in diagnostic or therapeutic procedures.

RGN Antibody (Center) Blocking Peptide - Protein Information

Name RGN

Synonyms SMP30

Function

Gluconolactonase with low activity towards other sugar lactones, including gulonolactone and galactonolactone. Can also hydrolyze diisopropyl phosphorofluoridate and phenylacetate (in vitro). Calcium-binding protein. Modulates Ca(2+) signaling, and Ca(2+)- dependent cellular processes and enzyme activities (By similarity).

Cellular Location

Cytoplasm.



RGN Antibody (Center) Blocking Peptide - Protocols

Provided below are standard protocols that you may find useful for product applications.

• Blocking Peptides

RGN Antibody (Center) Blocking Peptide - Images

RGN Antibody (Center) Blocking Peptide - Background

RGN is a highly conserved, calcium-binding protein, that is preferentially expressed in the liver and kidney. This protein may have an important role in calcium homeostasis. Studies in rat indicate that the protein may also play a role in aging, as it shows age-associated down-regulation.

RGN Antibody (Center) Blocking Peptide - References

Ishigami, A., Handa, S. Biosci. Biotechnol. Biochem. 67 (1), 158-160 (2003) Fujita, T., Shirasawa, T. Mech. Ageing Dev. 107 (3), 271-280 (1999) Fujita, T., Mandel, J.L. Biochim. Biophys. Acta 1263 (3), 249-252 (1995)