

Thromboxane B2 Blocking Peptide
Catalog # PBV10325b**Specification**

Thromboxane B2 Blocking Peptide - Product Information**Thromboxane B2 Blocking Peptide - Additional Information****Application & Usage**

The peptide is used for blocking the antibody activity of active Thromboxane B2. It usually blocks the antibody activity completely in Western blot analysis by incubating the peptide with equal volume of antibody for 30 minutes at 37°C

Target/Specificity
Thromboxane B2**Formulation**

50 µg (0.2 mg/ml) in phosphate buffered saline (PBS), pH 7.2, containing 0.1% BSA and 0.02% thimerosal.

Reconstitution & Storage
-20 °C**Background Descriptions****Precautions**

Thromboxane B2 Blocking Peptide is for research use only and not for use in diagnostic or therapeutic procedures.

Thromboxane B2 Blocking Peptide - Protein Information**Thromboxane B2 Blocking Peptide - Protocols**

Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)

Thromboxane B2 Blocking Peptide - Images

