

Human CellExp LIF, mouse recombinant protein

LIF, CDF, DIA, HILDA, MLPLI Catalog # PBV10893r

Specification

Human CellExp LIF, mouse recombinant protein - Product info

Primary Accession <u>P09056</u>

Calculated MW of 20.7 kDa with no tags.
The predicted N-terminus is Ser 24.
DTT-reduced Protein migrates as 33-45

kDa due to different glycosylation. KDa

Human CellExp LIF, mouse recombinant protein - Additional Info

Gene ID 16878
Gene Symbol LIF

Other Names

LIF, CDF, DIA, HILDA, MLPLI

Gene Source Mouse

Source HEK 293 cells
Assay&Purity SDS-PAGE; ≥92%

Assay2&Purity2 N/A;
Recombinant Yes

Target/Specificity

LIF

Application Notes

Centrifuge the vial prior to opening. Reconstitute in sterile deionized water to a concentration of $50 \mu g/ml$. Do not vortex. This solution can be stored at 2-8°C for up to 1 month. For extended storage, it is recommended to store at -20°C.

Format

Lyophilized powder

Storage

-20°C; Sterile filtered through a 0.22 micron filter. Lyophilized from $1 \times PBS$, pH 7.4. Generally 5-8% Mannitol or trehalose is added as a protectant before lyophilization.

Human CellExp LIF, mouse recombinant protein - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- Immunohistochemistry
- Immunofluorescence



- Immunoprecipitation
- Flow Cytomety
- Cell Culture

Human CellExp LIF, mouse recombinant protein - Images

Human CellExp LIF, mouse recombinant protein - Background

Leukemia inhibitory factor, or LIF, an interleukin 6 class cytokine, is a protein in cells that affects cell growth and development. Leukemia Inhibitory Factor has several functions such as cholinergic neuron differentiation, control of stem cell pluripotency, bone & fat metabolism, mitogenesis of factor dependent cell lines & promotion of megakaryocyte production in vivo. Removal of LIF pushes stem cells toward differentiation, but they retain their proliferative potential or pluripotency. Therefore LIF is used in mouse embryonic stem cell culture. It is necessary to maintain the stem cells in an undifferentiated state, however genetic manipulation of embryonic stem cells allows for LIF independent growth, notably overexpression of the gene Nanog. LIF is not required for culture of human embryonic stem cells.

Human CellExp LIF, mouse recombinant protein - References

Gough N.M.,et al.Nucleic Acids Res. 16:9857-9857(1988). Stahl J.,et al.J. Biol. Chem. 265:8833-8841(1990). Teuscher C.,et al.J. Immunol. 163:2262-2266(1999). Gearing D.P.,et al.EMBO J. 6:3995-4002(1987). Gough N.M.,et al.Blood Cells 14:431-442(1988).