

Theromacin (49 - 63)(C57-C59 disulfide bridge)
Synthetic Peptide
Catalog # SP2377a**Specification**

Theromacin (49 - 63)(C57-C59 disulfide bridge) - Product InformationPrimary Accession
Sequence[Q6T6C2](#)
NH2-QRLPNNKQCRCINAR -COOH**Theromacin (49 - 63)(C57-C59 disulfide bridge) - Additional Information****Other Names**

Theromacin, THMAC

Format

Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.

Storage

Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.

Precautions

This product is for research use only. Not for use in diagnostic or therapeutic procedures.

Theromacin (49 - 63)(C57-C59 disulfide bridge) - Protein Information**Name** THMAC**Function**

Has a bactericidal activity. Active against *M.luteus*. No activity toward *E.coli* and *F.oxysporum*.

Cellular Location

Secreted.

Tissue Location

Coelomic liquid (at protein level). Expressed in large fat cells in contact with coelomic cavities, in intestinal epithelia and at the epidermis level

Theromacin (49 - 63)(C57-C59 disulfide bridge) - Images